

Gas Chromatography Accessories

GC Ferrules

For Column Connections

GLScience's ferrules for gas chromatography by are manufactured from high purity graphite and graphite/Vespel, and are blister-packed to prevent contamination. Available in a variety of sizes, they are suitable for most GC and GC-MS.





GVL



GF



GVS

Ferrule for Agilent GC

Description	P/N	Applicable Column's I.D.	Remark	Qty.	Cat.No.
15% Graphite / 85% Vespel Ferrule	GVS-0.4	0.10~0.25mm	GC Injection Port, Detector	10pcs	3007-41140
	GVS-0.5	0.32mm		10pcs	3007-41150
	GVS-0.8	0.53mm		10pcs	3007-41180
15% Graphite / 85% Vespel Ferrule	GVL-0.4	0.10~0.25mm	Connect to GC/MS Interface	10pcs	3007-31144
	GVL-0.5	0.32mm		10pcs	3007-31145
	GVL-0.8	0.53mm		10pcs	3007-31148
Graphite Ferrule	GF-0.5	0.1~0.32mm	GC Injection Port, Detector	10pcs	3007-31305
	GF-0.8	0.53mm		10pcs	3007-31308
Super Graphite Ferrule	SGF-0.5	0.1~0.32mm	One-ring Ferrule	6pcs	3007-31405

Ferrule for Shimadzu GC


Description	P/N	Applicable Column's I.D.	Remark	Qty.	Cat.No.
15% Graphite / 85% Vespel Ferrule	GVL-0.4	0.10~0.25mm	GCMS Injection Port • MS Interface	10pcs	3007-31144
15% Graphite / 85% Vespel Ferrule	GVL-0.5	0.32mm		10pcs	3007-31145
15% Graphite / 85% Vespel Ferrule	GVL-0.8	0.53mm		10pcs	3007-31148
Graphite Ferrule	G-0.5	0.1~0.32mm	GC Injection, Detector	10pcs	3007-14005*
Graphite Ferrule	G-0.8	0.53mm		10pcs	3007-14008*
Super Graphite Ferrule	SG-0.5	0.1~0.32mm	Specific Design for 0.1~0.32mm I.D. Columns	6pcs	3007-14055

*It is not in Blister pack.

Vespel Ferrule Remover

Vespel ferrules stuck to the nut can be easily removed.

Description	Cat.No.
GL Remover	3001-12745



GL Remover